

Far Infrared Warming Products

Choose the best solution for your research



Kent Scientific's Far Infrared warming products provide a variety of options to safely warm your animals. From simple pad power regulation to automatic control based on temperature feedback, we have the solution to every research warming requirement.

Features & Benefits:

- **Deep Penetration**
Warms animal quickly and safely
- **Small Footprint**
Easily fits into your surgical setup
- **Homeothermic options**
Regulate animal temperature




Nylon fabric warming pad
available in 2 sizes

SOLE DISTRIBUTOR & AUTHORIZED REPRESENTATIVE


PELINDUNG MEDIC SDN BHD
(1274463-U)

9-2 Jalan USJ10/1G,
47620 Subang Jaya,
Selangor, Malaysia

+603-56338444
+6012-662 5371 / +6012-6709676
www.pelindungmedic.com
info@pelindungmedic.com









Kent Scientific
CORPORATION

Far Infrared Warming

Product Selection Guide



Warming Pad Control

- Small footprint; easily fits into your surgical set-up
- Large display; quickly set and identify target temperature
- Operates using AC power



RightTemp® Jr Homeothermic Control

- Automatic control measures and regulates your animal's temperature
- Includes temperature sensor
- Far infrared warming pad provided



PhysioSuite® RightTemp® Module with Advanced Features

- Temperature monitor and homeothermic control
- Modular physiological monitor system
- Optional Modules for ventilation, end-tidal CO₂, pulse oximetry and heart rate monitoring

	Warming Pad Control	RightTemp® Jr Homeothermic Control	PhysioSuite®
Homeothermic control	○	●	●
Pad temperature displayed	●	●	●
Body temperature displayed	○	●	●
Temperature sensor(s) included	●	●	●
15.2cm x 20.3cm Nylon fabric pad	●	●	○
20.3cm x 25.4cm Nylon fabric pad	●	●	●
Internal Memory	○	○	●
Data export capability	○	○	●
Real-time display	●	●	●
User custom menu	○	○	●
Adjustable temperature settings	○	●	●
AC powered	●	●	●
Electrophysiology compatible	○	○	●
Optional ventilator module	○	○	●
Optional CO ₂ module	○	○	●
Optional pulse oximetry module	○	○	●

SOLE DISTRIBUTOR & AUTHORIZED REPRESENTATIVE

PELINDUNG MEDIC SDN BHD
(1274463-U)

9-2 Jalan USJ10/1G,
47620 Subang Jaya,
Selangor, Malaysia

+603-56338444

+6012-662 5371 / +6012-6709676

www.pelindungmedic.com

info@pelindungmedic.com

CGMHG
Good Distribution
Practice for
Medical Devices

Medical Device
AUTHORITY
MALAYSIA



Kent Scientific
CORPORATION